

ReMove: A Configuration Tool for Robot-Assisted Experiments in Rehabilitation

Qiang Xu^{*3}, Zhi Chen^{†1, 2}, Jiayi Tang³, Cheng Zhuang³, Teck Hong Koh³, and Chuanxin M. Niu^{‡1, 4}

¹ Department of Rehabilitation Medicine, Ruijin Hospital, School of Medicine, Shanghai Jiao Tong University, Shanghai, China. ² School of Medicine, Shanghai Jiao Tong University, Shanghai, China. ³ Shanghai Fourier Intelligence Co., Ltd. ⁴ Laboratory of Neurorehabilitation Engineering, School of Biomedical Engineering, Shanghai Jiao Tong University, Shanghai, China.

DOI: [10.21105/joss.04071](https://doi.org/10.21105/joss.04071)

Software

- [Review](#)
- [Repository](#)
- [Archive](#)

Editor: [Adi Singh](#)

Reviewers:

- [@AustinTSchaffer](#)
- [@sneakers-the-rat](#)

Submitted: 16 December 2021

Published: 17 January 2022

License

Authors of papers retain copyright and release the work under a Creative Commons Attribution 4.0 International License ([CC BY 4.0](#)).

Summary

Resistance training has shown efficacy in post-stroke rehabilitation. To understand which aspects of resistance (i.e. magnitude, type) play an important role in movement rehabilitation, haptic robots are widely used for simulation of resistance due to easy tweaking of mechanical parameters. Although adjustment of robot-delivered resistance can be achieved by direct modification on the source code, it would be much more efficient if the configuration can be integrated into a standalone toolkit. Therefore we have developed ReMove, a programmable application for customization of robot-delivered resistance and robot-guided movements on rehabilitation robots. ReMove provides a plain-text interface that specifies the robot behavior during each trial of movement, and it renders both the graphical and mechanical specs during the movement. By the time of submission, ReMove had been supporting 2 clinical studies. Pilot data supported the feasibility and utility of ReMove in movements against light and heavy resistance.

Statement of need

Previous studies have shown the therapeutic effect of resistance training, including the improvement of gait speed (Mehta et al. 2012), muscle strength (U Flansbjerg, Lexell, and Brogårdh 2012), perceived participation (U Flansbjerg et al. 2008), etc. It becomes interesting to test which aspects of mechanical resistance could elicit physiological responses that are beneficial. This line of studies requires easy tweaking of resistance, including the change of type, magnitude, and timing of resistance. A major challenge for this, however, is the generation of different types of resistance, e.g. it is cumbersome to prepare viscosity in real-world setups. The second challenge arises in flexibly changing the parameters using real-world objects, e.g. alternating between heavy and light resistance across trials.

Simulation of resistance using haptic robots are advantageous for learning the sensorimotor control in humans. The properties of resistance can be changed in the source code on demand. The disadvantage is also obvious that most clinical researchers lack the knowledge of programming. Although the barrier may be lowered for robot programming by use of graphical environments (LabVIEW, Simulink) or tools for experiment customization (E-Prime,

^{*}co-first author

[†]co-first author

[‡]corresponding author: minos.niu@gmail.com

37 PsychoPy), no existing software caters to the need of simulation and coordination of resistance
38 in rehabilitation studies.

39 Here we introduce ReMove, a programmable application for the customization of robot-
40 delivered resistance and robot-guided movements on haptic robots. The main value of ReMove
41 is to provide an integrated research environment that allowed both generation and manipula-
42 tion of resistance during human upper-limb movements. Using ReMove, it would significantly
43 reduce the time to prepare an experiment that needs trial-by-trial modification on resistance.

44 Software overview

45 The architecture of ReMove is shown in Figure 1. Notice that the experimental setups are
46 specified through the Parameter Interfaces, and most of the logics are programmed in the
47 module of device synchronization. Modules communicate asynchronously for better perfor-
48 mance. The current version of ReMove is adapted to M2PRO, a 2D haptic robot licensed for
49 clinical use in China (Fourier Intelligence inc., Shanghai, China). The support will expand to
50 other available makes and models including Phantom series, HapticMASTER, etc.





Figure 1: The architecture of ReMove.

51 ReMove allows researchers to customize the following parameters using plain-text (.csv format
52 in the existing version):

- 53 1. The location of the starting point.
- 54 2. The location of the target.
- 55 3. The type of resistance.
- 56 4. The magnitude of resistance.

57 During each trial of movement, ReMove renders the graphical details corresponding to the .csv
58 files. Both a subject view (Figure 2A, minimal display of experimental status) an inspector
59 view (Figure 2B, rich information about experimental status) are rendered.

A)



B)




Figure 2: The main UIs of ReMove. A) the subject view, displaying only task-related elements to the participant. B) the inspector view, with additional information about test progress, device status, etc.

Experimental example

61 We accomplished pilot experiments with a volunteer. The volunteer was asked to move from a
62 starting point (close to his chest) to 7 different targets (25cm away) under 2 levels of inertial
63 load. Parameter setting was accomplished by a therapist naïve of computer programming.
64 The parameter setting took about 5 minutes. Figure 3 shows that the volunteer performed
65 straight movements. Notice that the peak-velocity decreased (Figure 3B) due to the increased
66 magnitude of resistance.




Figure 3: The Pilot data of ReMove. A) The actual scene of one clinical study supported by ReMove. B) This panel displays the trajectories and movement velocity profiles in light condition ($10\text{Ns}^2/\text{m}$, $30\text{Ns}/\text{m}$). C) This panel displays the trajectories and movement velocity profiles in heavy condition ($50\text{Ns}^2/\text{m}$, $30\text{Ns}/\text{m}$).

Acknowledgements

The authors thank Yongjun Qiao for help with the design and test of ReMove. This work is funded by the General Project of National Natural Science Foundation of China (81971722); Shanghai Science and Technology Commission Project (19511105600); Shanghai Municipal Health Commission Project (2019SY004).

References